



Application Note AN N308

Analysis of Pet Food with FT-NIR Spectroscopy

Pets are today an important part of people's lives. Therefore consumers trust in brands which guarantee the best quality for their individual animal, since pets of different breed, size and age often need a special nutrient profile in their diet. Pet food producers have to provide a vast amount of different products, produced at an acceptable cost.

From raw materials to finished pet food

The biggest challenge in pet food production is to keep track of the different recipes depending on the availability of the expensive raw materials. Depending on origin and season, they often show huge batch-to-batch variability; nevertheless, feeding the least cost formulation programs with valid compositional data is essential and a fast and accurate analysis, providing results in seconds can help keeping track of the production process.

Near infrared spectroscopy is used in agriculture for many decades for the analysis of main constituents like moisture, fat and protein. It is fast and safe, since it requires little to no sample preparation, nor hazardous chemicals or gases. This makes a quick test of the incoming raw material using FT-NIR available even to production staff, before processing the materials and also allows a fast way of testing the final material before packaging.

Easy sampling and measurement

Solid samples like raw materials and finished pet food are filled in cups with a quartz glass bottom, placed on the integrating sphere of the MPA II or TANGO-R spectrometer and measured from below. This way, the product has an even surface and for highly reproducible results, it is rotated during the measurement. Liquid samples like oils and fats are simply filled in 8mm glass vials and are measured in the sample compartment of the MPA II or TANGO-T.



Analyzing rabbit food in diffuse reflection on the integrating sphere of the MPA II FT-NIR spectrometer.

Analyzing fat quality


The fat content in pet food mainly derives from the protein source (e.g. the fat in poultry) but also from added fats e.g. lard, tallow, vegetable or fish oil. Quality pet foods will use good sources of fat and ensure they provide an optimized level of omega-3 oils. Here FT-NIR spectroscopy can help to determine the amount of EPA, DHA and total omega-3 in fish oils. Moreover, the fat quality in general in terms of trans fatty acids (TFA), free fatty acids (FFA) and fatty acid composition can be analyzed.

Online analysis for real-time monitoring

Even real-time online measurements can be carried out using FT-NIR process spectrometers, ensuring highest product quality at all times and minimized rework to maximize profits. Depending on the output of the factory, a return on investment on NIR technology can be realized within a year.

Ready to use calibration packages

A comprehensive set of universal calibrations for the analysis of raw materials and finished products in the pet food industry as well as for various edible oils and fats is available on request. These calibrations enable a quick and cost-effective start to your FT-NIR analysis, in the lab or online.



Calibration graphs for protein and fat with FT-NIR spectroscopy.

Measurable Parameters in Petfood:

- Moisture
- Fat
- Protein
- Fiber
- Ash
- Starch

FT-NIR Spectrometers: Bruker Optics offers various FT-NIR spectrometer models for lab, at-line and on-line applications:

TANGO



FT-NIR analyzer for routine use in the lab.

MPA II



Multi Purpose Analyzer for maximum flexibility.

MATRIX-I



At-line analysis with optional NEMA4/IP66 protection.

MATRIX-F



Process monitoring with probes and measurement heads.

● Bruker Optik GmbH

Ettlingen · Deutschland
Phone +49 (7243) 504-2000
Fax +49 (7243) 504-2050
info.bopt.de@bruker.com

www.bruker.com/optics

Bruker Optics is continually improving its products and reserves the right to change specifications without notice.
© 2018 Bruker Optics BOPT-4001224-01

Bruker Optics Inc.

Billerica, MA · USA
Phone +1 (978) 439-9899
Fax +1 (978) 663-9177
info.bopt.us@bruker.com

Bruker Shanghai Ltd.

Shanghai · China
Phone +86 21 51720-890
Fax +86 21 51720-899
info.bopt.cn@bruker.com